

An Open Pediatric Brain Tumor Atlas

This manuscript ([permalink](#)) was automatically generated from [AlexsLemonade/OpenPBTA-manuscript@ac9b7cd](#) on September 16, 2019.

Authors

- **John Doe**

 [XXXX-XXXX-XXXX-XXXX](#) ·  [johndoe](#) ·  [johndoe](#)

Department of Something, University of Whatever · Funded by Grant XXXXXXXX

- **Jane Roe**

 [XXXX-XXXX-XXXX-XXXX](#) ·  [janeroe](#)

Department of Something, University of Whatever; Department of Whatever, University of Something

Abstract

Introduction

Introduction will go here.

Materials and Methods

Biospecimen collection

The Pediatric Brain Tumor Atlas specimens are comprised of samples from Children's Brain Tumor Tissue Consortium (CBTTC) and the Pediatric Pacific Neuro-oncology Consortium (PNOC).

Children's Brain Tumor Tissue Consortium (CBTTC)

The CBTTC [1] is a collaborative, multi-institutional (16 institutions worldwide) research program dedicated to the study of childhood brain tumors. All CBTTC data can be download from the Gabriella Miller Kids First Data Resource Center (KF-DRC), [2]. The deidentified patient's blood and tumor tissue were prospectively collected by the consortium from patients enrolled within the CBTTC.

The cell lines were generated by the CBTTC from either fresh tumor tissue obtained directly from surgery performed at Children's Hospital of Philadelphia (CHOP) or from prospectively collected tumor specimens stored in Recover Cell Culture Freezing media (cat# 12648010, Gibco). The tissue was dissociated using enzymatic method with papain as described [3]. Briefly, tissue was washed with HBSS (cat# 14175095, Gibco), minced and incubated with activated papain solution (cat# LS003124, SciQuest) for up to 45 minutes. The papain was inactivated using ovomucoid solution (cat# 542000, SciQuest), tissue was briefly treated with DNase (cat# 10104159001, Sigma) and passed through the 100µm cell strainer (cat# 542000, Greiner Bio-One). Two cell culture conditions were initiated based on the number of cells available. For cultures utilizing the fetal bovine serum (FBS), a minimum density of 3×10^5 cells/ml were plated in DMEM/F-12 medium (cat# D8062, Sigma) supplemented with 20% FBS (cat# SH30910.03, Hyclone), 1% GlutaMAX (cat# 35050061, Gibco), Penicillin/Streptomycin-Amphotericin B Mixture (cat# 17-745E, Lonza) and 0.2% Normocin (cat# ant-nr-2, Invivogen). For the serum-free media conditions cells were plated at minimum density of 1×10^6 cells/ml in DMEM/F12 media supplemented with 1% GlutaMAX, 1x B-27 supplement minus vitamin A (cat# 12587-010, Gibco), 1x N-2 supplement (cat# 17502001, Gibco), 20 ng/ml epidermal growth factor (cat# PHG0311L, Gibco), 20 ng/ml basic fibroblast growth factor (cat# 100-18B, PeproTech), 2.5µg/ml heparin (cat# H3149, Sigma), Penicillin/Streptomycin-Amphotericin B Mixture and 0.2% Normocin.

Pacific Pediatric Neuro Oncology Consortium (PNOC)

The Pacific Pediatric Neuro-Oncology Consortium (PNOC) is an international consortium dedicated to bringing new therapies to children and young adults with brain tumors. PNOC collected blood and tumor biospecimens from newly-diagnosed DIPG patients as part of the clinical trial [PNOC003/NCT02274987](#) [4].

Nucleic acids extraction and library preparation

PNOC samples

The Translational Genomic Research Institute (TGEN; Phoenix, AZ) performed DNA and RNA extractions on tumor biopsies using a DNA/RNA AllPrep Kit (Qiagen, #80204). All RNA used for library prep had a minimum RIN of 7 but no QC thresholds were implemented for the DNA. For library

preparation, 500ng of nucleic acids were used as input for RNA-Seq and WXS. The RNA prep was performed using the TruSeq RNA Sample Prep Kit (Illumina, #FC-122-1001) and the exome prep was performed using KAPA Library Preparation Kit (Kapa Biosystems, #KK8201) using Agilent's SureSelect Human All Exon V5 backbone with custom probes. These probes include CGH probes that target 44,000 evenly spaced genomic loci to assess copy number changes, as well as probes that tile across tumor suppressor genes and genes involved in common cancer translocations. All extractions and library preparations were performed according to manufacturer's instructions.

CBTTC samples

Blood, tissue and cell line DNA/RNA extraction was performed at Biorepository Core (BioRC) at CHOP. Briefly, 10-20 mg frozen tissue, 0.4-1ml of blood or 2×10^6 cells pellet was used for extractions. Tissues were lysed using a Qiagen TissueLyser II (Qiagen) with 2×30 sec at 18Hz settings using 5 mm steel beads (cat# 69989, Qiagen). Both tissue and cell pellets processes included a CHCl₃ extraction and were run on the QiaCube automated platform (Qiagen) using the AllPrep DNA/RNA/miRNA Universal kit (cat# 80224, Qiagen). Blood was thawed and treated with RNase A (cat#, 19101, Qiagen); 0.4-1ml was processed using the Qiagen QIAasympy automated platform (Qiagen) using the QIAasympy DSP DNA Midi Kit (cat# 937255, Qiagen). DNA and RNA quantity and quality was assessed by PerkinElmer DropletQuant UV-VIS spectrophotometer (PerkinElmer) and an Agilent 4200 TapeStation (Agilent, USA) for RINe and DINe (RNA Integrity Number equivalent and DNA Integrity Number equivalent respectively). Library preparation and sequencing was performed by the NantHealth sequencing center. Briefly, DNA sequencing libraries were prepared for tumor and matched-normal DNA using the KAPA Hyper prep kit (cat# KK8541, Roche); tumor RNA-Seq libraries were prepared using KAPA Stranded RNA-Seq with RiboErase kit (cat# KK8484, Roche). Whole genome sequencing (WGS) was performed at an average depth of coverage of 60X for tumor samples and 30X for germline. RNA samples were sequenced to an average of 200M reads. All samples were sequenced on the Illumina HiSeq platform (X/400) (Illumina) with 2×150 bp read length. For the cell line sequencing, samples labelled with "CL-adh" correspond to the adherent FBS cell lines and those labelled "CL-susp" are the serum-free lines.

Data generation

NantHealth Sequencing Center (Culver City, CA) performed whole genome sequencing (WGS) on all paired tumor (~60X) and constitutive (~30X) DNA samples. WGS libraries were 2×150 bp and sequenced on an Illumina X/400. NantHealth Sequencing Center performed ribosomal-depleted whole transcriptome stranded RNA-Seq to an average depth of 100M reads for CBTTC tumor samples. The Translational Genomic Research Institute (TGEN; Phoenix, AZ) performed paired tumor (~200X) and constitutive whole exome sequencing (WXS) and poly-A selected RNA-Seq (~200M reads) for PNOc tumor samples. PNOc WXS and RNA-Seq libraries 2×100 bp and sequenced on an Illumina HiSeq 2500.

DNA WGS Alignment

We used BWA-MEM [5] v0.7.17 for alignment of paired-end DNA-seq reads. The alignment reference that we used was Homo Sapiens Human Genome (hg) version 38, patch release 12, fasta file obtained from UCSC [6]. Alignments were further processed using following the Broad Institute's Best Practices [7] for processing BAMs in preparation for variant discovery. Duplicates were marked using Samblaster[8] v0.1.24, BAMs merged and sorted using Sambamba [9] v0.6.3. Lastly, resultant BAMs were processed using Broad's Genome Analysis Tool Kit (GATK) [10] v4.0.3.0, BaseRecalibrator submodule.

Germ Line Single Nucleotide Variant Calling

Somatic Single Nucleotide Variant Calling

SNV and INDEL calling

We used Strelka2 [[11](#)] v2.9.3 and Mutect2 from GATK v4.1.1.0. Strelka2 was run using default parameters on human genome reference hg38, canonical chromosomes only (chr1-22, X,Y,M), as recommended by the author. Mutect2 was run following Broad best practices outlined from their Workflow Description Language (WDL) [[12](#)].

VCF annotation and MAF creation

We filtered outputs from both callers on the “PASS” filter, and annotated using The ENSEMBL Variant Effect Predictor [[13](#)], reference release 93, and created MAFs using MSKCC’s vcf2maf [[14](#)] v1.6.16.

Somatic Copy Number Variant Calling

We used Control-FREEC [[15](#), [16](#)] v8.7 and CNVkit [[17](#)] v0.9.3 for copy number variant calls. CNVkit was run using default parameters on human genome reference hg38 and using the batch command for tumor-normal pairs rather than a panel of normals.

Somatic Structural Variant Calling

We used Manta SV [[18](#)] v1.4.0 for structural variant (SV) calls. Manta SV calling was also limited to regions used in Strelka2. We also ran LUMPY SV [[19](#)] v0.2.13 in express mode using default parameters. The hg38 reference used was also limited to canonical chromosome regions.

Gene Expression Abundance Estimation

We used STAR [[20](#)] v2.6.1d to align paired-end RNA-seq reads. This output was used for all subsequent RNA analysis. The reference we used was that of ENSEMBL’s GENCODE 27 [[21](#)], “Comprehensive gene annotation.” We used RSEM [[22](#)] v1.3.1 for transcript- and gene-level quantification. We also added a second method of quantification using kallisto [[23](#)] v0.43.1. This method differs in that it uses pseudoalignments using fastq reads directly to the aforementioned GENCODE 27 reference.

RNA Fusion Calling and Prioritization

Gene fusion detection

We set up [Arriba v1.1.0](#) and STAR-Fusion 1.5.0 [[24](#)] fusion detection tools using CWL on CAVATICA. For both these tools we used aligned BAM and chimeric SAM files from STAR as inputs and GRCh38_gencode_v27 GTF for gene annotation. We ran STAR-Fusion with default parameters and annotated all fusion calls with GRCh38_v27_CTAT_lib_Feb092018.plugin-play.tar.gz provided in the STAR-fusion release. For Arriba, we used a blacklist file (blacklist_hg38_GRCh38_2018-11-04.tsv.gz) from the Arriba release tarballs to remove recurrent fusion artifacts and transcripts present in healthy tissue. We also provided Arriba with strandedness information or set it to auto-detection for polyA samples.

Fusion prioritization

We built a [fusion prioritization pipeline](#) to filter and annotate fusions. We considered all inframe and frameshift fusion calls with 1 or more junction reads and fused genes expressed with TPM greater than one to be true calls. If a fusion call had large number of spanning fragment reads compared to

junction reads (spanning fragment minus junction read greater than ten) or if either 5' or 3' genes fused to more than five different genes we removed these calls as a potential false positive. We also removed fusions if the 5' or 3' ends were the same gene, and these were tagged as non-canonical splicing or duplication. We used a list of curated fusion calls for each histology to capture each occurrence of the fusion as a putative driver fusion. We prioritized a union of fusion calls as true calls if the fused genes were detected by both callers, the same fusion was recurrent in histology (>2 samples) or the fusion was specific to the broad histology. We annotated putative driver fusions and prioritized fusions lists with kinases, oncogenic, tumor suppressor, transcription factor, fused genes and known TCGA fusions from curated [datasheets](#). We also added chimerDB [\[25\]](#) annotations to both driver and prioritized fusion list.

Clinical Data Harmonization

Prediction of participants' genetic sex

The clinical metadata provided included a reported gender. We used DNA data, in concert with the reported gender, to predict participant genetic sex so that we could identify sexually dimorphic outcomes. This analysis could also reveal samples that may have been contaminated in certain circumstances. We used the idxstats utility from SAMTOOLS [\[26\]](#) to calculate read lengths, the number of mapped reads, and the corresponding chromosomal location for reads to the X and Y chromosomes. We used the fraction of total normalized X and Y chromosome reads that were attributed to the Y chromosome as a summary statistic. We reviewed this statistic in the context of reported gender and determined that a threshold of less than 0.2 clearly delineated female samples. Fractions greater than 0.4 were predicted to be males. Samples with values in the range [0.2, 0.4] were marked as unknown. We ran this analysis through [CWL](#) on Cavatica. Resulting calls were added to the clinical metadata as `germline_sex_estimate`.

Results

Results section stub.

Conclusions

Stub in conclusions section

References

1. Home

Children's Brain Tumor Tissue Consortium

<https://cbttc.org/>

2. Working Together to Put Kids First <https://kidsfirstdrc.org/>

3. Pediatric High Grade Glioma Resources From the Children's Brain Tumor Tissue Consortium (CBTTC) and Pediatric Brain Tumor Atlas (PBTA)

Heba Ijaz, Mateusz Koptyra, Krutika S. Gaonkar, Jo Lynne Rokita, Valerie P. Baubet, Lamiya Tauhid, Yankun Zhu, Miguel Brown, Gonzalo Lopez, Bo Zhang, ...

Cold Spring Harbor Laboratory (2019-05-31) <https://doi.org/gf66qt>

DOI: [10.1101/656587](https://doi.org/10.1101/656587)

4. A pilot precision medicine trial for children with diffuse intrinsic pontine glioma—PNOC003: A report from the Pacific Pediatric Neuro-Oncology Consortium

Sabine Mueller, Payal Jain, Winnie S. Liang, Lindsay Kilburn, Cassie Kline, Nalin Gupta, Eshini Panditharatna, Suresh N. Magge, Bo Zhang, Yuankun Zhu, ... Adam C. Resnick

International Journal of Cancer (2019-04-03) <https://doi.org/gf6pfb>

DOI: [10.1002/ijc.32258](https://doi.org/10.1002/ijc.32258) · PMID: [30861105](https://pubmed.ncbi.nlm.nih.gov/30861105/)

5. Aligning sequence reads, clone sequences and assembly contigs with BWA-MEM

Heng Li

arXiv (2013-03-16) <https://arxiv.org/abs/1303.3997v2>

6. Index of /goldenPath/hg38/bigZips <http://hgdownload.soe.ucsc.edu/goldenPath/hg38/bigZips/>

7. GATK | BP Doc # | <https://software.broadinstitute.org/gatk/best-practices/workflow?id>

8. SAMBLASTER: fast duplicate marking and structural variant read extraction

G. G. Faust, I. M. Hall

Bioinformatics (2014-05-07) <https://doi.org/f6kft3>

DOI: [10.1093/bioinformatics/btu314](https://doi.org/10.1093/bioinformatics/btu314) · PMID: [24812344](https://pubmed.ncbi.nlm.nih.gov/24812344/) · PMCID: [PMC4147885](https://pubmed.ncbi.nlm.nih.gov/PMC4147885/)

9. Sambamba: fast processing of NGS alignment formats

Artem Tarasov, Albert J. Vilella, Edwin Cuppen, Isaac J. Nijman, Pjotr Prins

Bioinformatics (2015-02-19) <https://doi.org/gfzsfw>

DOI: [10.1093/bioinformatics/btv098](https://doi.org/10.1093/bioinformatics/btv098) · PMID: [25697820](https://pubmed.ncbi.nlm.nih.gov/25697820/) · PMCID: [PMC4765878](https://pubmed.ncbi.nlm.nih.gov/PMC4765878/)

10. GATK | Home <https://software.broadinstitute.org/gatk/>

11. Strelka2: fast and accurate calling of germline and somatic variants

Sangtae Kim, Konrad Scheffler, Aaron L. Halpern, Mitchell A. Bekritsky, Eunho Noh, Morten Källberg, Xiaoyu Chen, Yeonbin Kim, Doruk Beyter, Peter Krusche, Christopher T. Saunders

Nature Methods (2018-07-16) <https://doi.org/gdwrp4>

DOI: [10.1038/s41592-018-0051-x](https://doi.org/10.1038/s41592-018-0051-x) · PMID: [30013048](https://pubmed.ncbi.nlm.nih.gov/30013048/)

12. Official code repository for GATK versions 4 and up: broadinstitute/gatk

Broad Institute

(2019-09-16) <https://github.com/broadinstitute/gatk>

13. The Ensembl Variant Effect Predictor

William McLaren, Laurent Gil, Sarah E. Hunt, Harpreet Singh Riat, Graham R. S. Ritchie, Anja Thormann, Paul Flicek, Fiona Cunningham

Genome Biology (2016-06-06) <https://doi.org/gdz75c>

DOI: [10.1186/s13059-016-0974-4](https://doi.org/10.1186/s13059-016-0974-4) · PMID: [27268795](https://pubmed.ncbi.nlm.nih.gov/27268795/) · PMCID: [PMC4893825](https://pubmed.ncbi.nlm.nih.gov/PMC4893825/)

14. Convert a VCF into a MAF, where each variant is annotated to only one of all possible gene isoforms: mskcc/vcf2maf

Memorial Sloan Kettering

(2019-09-10) <https://github.com/mskcc/vcf2maf>

15. Control-FREEC: a tool for assessing copy number and allelic content using next-generation sequencing data

Valentina Boeva, Tatiana Popova, Kevin Bleakley, Pierre Chiche, Julie Cappel, Gudrun Schleiermacher, Isabelle Janoueix-Lerosey, Olivier Delattre, Emmanuel Barillot

Bioinformatics (2011-12-06) <https://doi.org/ckt4vz>

DOI: [10.1093/bioinformatics/btr670](https://doi.org/10.1093/bioinformatics/btr670) · PMID: [22155870](https://pubmed.ncbi.nlm.nih.gov/22155870/) · PMCID: [PMC3268243](https://pubmed.ncbi.nlm.nih.gov/PMC3268243/)

16. Control-free calling of copy number alterations in deep-sequencing data using GC-content normalization

Valentina Boeva, Andrei Zinovyev, Kevin Bleakley, Jean-Philippe Vert, Isabelle Janoueix-Lerosey, Olivier Delattre, Emmanuel Barillot

Bioinformatics (2010-11-15) <https://doi.org/c6bcps>

DOI: [10.1093/bioinformatics/btq635](https://doi.org/10.1093/bioinformatics/btq635) · PMID: [21081509](https://pubmed.ncbi.nlm.nih.gov/21081509/) · PMCID: [PMC3018818](https://pubmed.ncbi.nlm.nih.gov/PMC3018818/)

17. CNVkit: Genome-Wide Copy Number Detection and Visualization from Targeted DNA Sequencing

Eric Talevich, A. Hunter Shain, Thomas Botton, Boris C. Bastian

PLOS Computational Biology (2016-04-21) <https://doi.org/c9pd>

DOI: [10.1371/journal.pcbi.1004873](https://doi.org/10.1371/journal.pcbi.1004873) · PMID: [27100738](https://pubmed.ncbi.nlm.nih.gov/27100738/) · PMCID: [PMC4839673](https://pubmed.ncbi.nlm.nih.gov/PMC4839673/)

18. Manta: rapid detection of structural variants and indels for germline and cancer sequencing applications

Xiaoyu Chen, Ole Schulz-Trieglaff, Richard Shaw, Bret Barnes, Felix Schlesinger, Morten Källberg, Anthony J. Cox, Semyon Kruglyak, Christopher T. Saunders

Bioinformatics (2015-12-08) <https://doi.org/gf3ggb>

DOI: [10.1093/bioinformatics/btv710](https://doi.org/10.1093/bioinformatics/btv710) · PMID: [26647377](https://pubmed.ncbi.nlm.nih.gov/26647377/)

19. LUMPY: a probabilistic framework for structural variant discovery

Ryan M Layer, Colby Chiang, Aaron R Quinlan, Ira M Hall

Genome Biology (2014) <https://doi.org/gf3ggc>

DOI: [10.1186/gb-2014-15-6-r84](https://doi.org/10.1186/gb-2014-15-6-r84) · PMID: [24970577](https://pubmed.ncbi.nlm.nih.gov/24970577/) · PMCID: [PMC4197822](https://pubmed.ncbi.nlm.nih.gov/PMC4197822/)

20. STAR: ultrafast universal RNA-seq aligner

Alexander Dobin, Carrie A. Davis, Felix Schlesinger, Jorg Drenkow, Chris Zaleski, Sonali Jha, Philippe Batut, Mark Chaisson, Thomas R. Gingeras

Bioinformatics (2012-10-25) <https://doi.org/f4h523>

DOI: [10.1093/bioinformatics/bts635](https://doi.org/10.1093/bioinformatics/bts635) · PMID: [23104886](https://pubmed.ncbi.nlm.nih.gov/23104886/) · PMCID: [PMC3530905](https://pubmed.ncbi.nlm.nih.gov/PMC3530905/)

21. GENCODE - Human Release 27 https://www.gencodegenes.org/human/release_27.html

22. RSEM: accurate transcript quantification from RNA-Seq data with or without a reference genome

Bo Li, Colin N Dewey

BMC Bioinformatics (2011-08-04) <https://doi.org/cwg8n5>

DOI: [10.1186/1471-2105-12-323](https://doi.org/10.1186/1471-2105-12-323) · PMID: [21816040](https://pubmed.ncbi.nlm.nih.gov/21816040/) · PMCID: [PMC3163565](https://pubmed.ncbi.nlm.nih.gov/PMC3163565/)

23. Near-optimal probabilistic RNA-seq quantification

Nicolas L Bray, Harold Pimentel, Páll Melsted, Lior Pachter

Nature Biotechnology (2016-04-04) <https://doi.org/f8nvsp>

DOI: [10.1038/nbt.3519](https://doi.org/10.1038/nbt.3519) · PMID: [27043002](https://pubmed.ncbi.nlm.nih.gov/27043002/)

24. STAR-Fusion: Fast and Accurate Fusion Transcript Detection from RNA-Seq

Brian J. Haas, Alex Dobin, Nicolas Stransky, Bo Li, Xiao Yang, Timothy Tickle, Asma Bankapur, Carrie Ganote, Thomas G. Doak, Nathalie Pochet, ... Aviv Regev

Cold Spring Harbor Laboratory (2017-03-24) <https://doi.org/gf5pc5>

DOI: [10.1101/120295](https://doi.org/10.1101/120295)

25. OUP accepted manuscript

Nucleic Acids Research

(2016) <https://doi.org/gf6bx9>

DOI: [10.1093/nar/gkw1083](https://doi.org/10.1093/nar/gkw1083) · PMID: [27899563](https://pubmed.ncbi.nlm.nih.gov/27899563/) · PMCID: [PMC5210563](https://pubmed.ncbi.nlm.nih.gov/PMC5210563/)

26. The Sequence Alignment/Map format and SAMtools.

Heng Li, Bob Handsaker, Alec Wysoker, Tim Fennell, Jue Ruan, Nils Homer, Gabor Marth, Goncalo Abecasis, Richard Durbin,

Bioinformatics (Oxford, England) (2009-06-08) <https://www.ncbi.nlm.nih.gov/pubmed/19505943>

DOI: [10.1093/bioinformatics/btp352](https://doi.org/10.1093/bioinformatics/btp352) · PMID: [19505943](https://pubmed.ncbi.nlm.nih.gov/19505943/) · PMCID: [PMC2723002](https://pubmed.ncbi.nlm.nih.gov/PMC2723002/)